












Product Family
Titmus V4 Vision Screener

Select a Slide Set to form your complete instrument.
Refer to SKUs and Downloads for more information.



Product Numbers & Ordering Information

Product Numbers	Details
 <p>19108</p>	<p>Titmus V4 TNO Occupational Model Vision Screener with tests typically used by Occupational Health Doctors/Nurses for screening employees for various job functions, DOT testing. Tests include Visual Acuity (using Landolt rings per ANSI/ISO Stds), Binocular Vision, Depth / Color Perception, Muscle Balance (horizontal/vertical) and peripheral vision.</p>
 <p>19086</p>	<p>Titmus V4 Occupational Model Vision Screener with tests typically used by Occupational Health Doctors/Nurses for screening employees for various job functions, DOT testing. Tests include Visual Acuity (using Landolt rings/unbroken ring), Binocular Vision, Depth / Color Perception, Muscle Balance (horizontal/vertical) and peripheral vision.</p>
 <p>19089</p>	<p>Titmus V4 Occupational (Letters) Vision Screener with tests typically used by Occupational Health Doctors/Nurses for screening employees for various job functions, DOT testing. Tests include Visual Acuity (using Snellen letters), Binocular Vision, Depth / Color Perception, Muscle Balance (horizontal/vertical) and peripheral vision.</p>
 <p>19099</p>	<p>Titmus V4 Professional Model Vision Screener with tests typically used by Family Physicians, General Practitioners for general vision screening of adults and children. Tests include Visual Acuity (using letters/objects for adults/children), Plus Lens (children), Binocular Vision (children), Color Perception, Muscle Balance (horizontal/vertical) and peripheral vision.</p>
 <p>19090</p>	<p>Titmus V4 Aeromedical Model Vision Screener with tests typically used by AMEs (Aeromedical Examiner), General Practitioners for testing pilots' vision. Tests include Visual Acuity (using letters), Acuity at 32", Color Perception, Muscle Balance (horizontal/vertical) and peripheral vision.</p>
 <p>19097</p>	<p>Titmus V4 Pediatric Model Vision Screener with tests typically used by Pediatricians, General Practitioners for general vision screening of children. Tests per AAP recommendations include Visual Acuity (Michigan Preschool Test, Snellen E, Letters), Plus Lens (children), Binocular Vision, Color Perception, Muscle Balance and peripheral vision.</p>

	19100	<p>Titmus V4 School Model Vision Screener with tests typically used by General Practitioners, Pediatricians for general vision screening of school children. Tests include Visual Acuity (using letters/Snellen E), Plus Lens (children), Binocular Vision, Color Perception, Muscle Balance (horizontal/vertical) and peripheral vision.</p>
	19098	<p>Titmus V4 Preschool Model Vision Screener with tests typically used by General Practitioners, Pediatricians for general vision screening of preschool children. Tests include Visual Acuity (Efron Preschool Test, Snellen E, Letters), Plus Lens (children), Binocular Vision, Muscle Balance and peripheral vision.</p>
	19091	<p>Titmus V4 Driver Ed Model Vision Screener with tests typically used by General Practitioners for general vision screening of drivers. Tests include Visual Acuity (Letters, Numbers), Color Identification (traffic signal), Sign Recognition, Depth Perception, Muscle Balance (horizontal/vertical) and peripheral vision.</p>

Packaging Information

Packed 1 each in a cardboard carton.
 Packaged dimensions =
 Packaged weight = 22 lbs (10kg)
 Packaged dimensions=W x L x H 15.25" x 24" x 12.25" 38cm x 60cm x 31cm

Hazardous Materials

None

Overview

Key Features

The Titmus V4 vision screening instrument features a compact lightweight foldable design for ease of handling, membrane panel for operation, ergonomic design for patient comfort, standardized test illumination and an integrated handle for ease of transport.

- Test Distances:
- Far (20 feet/6m)
 - Near (14 inches/35cm)
 - Intermediate Distance (Optional)

Recommended Industries/Use

- Medical

Regulations

- EN 60601-2 - Meets CE Standards: Europe: EN 55011 Group 1, Class A / EN 60601-2 IEC 601-1-2
- UL 2601-1 - Meets UL Standards: USA: UL 2601-1, 1st Ed. Category 355
- CSA-C22.2 - Meets CSA Standards: Canada: CSA-C22.2, No. 601.1-M90, Category 245, Class 1 Device

Historical Brand

Titmus
 V4 replaces the following earlier models: Titmus i400, Titmus i500, Titmus 2A, and Titmus IIR

Warranty Information

Product has a warranty for a period of two (2) years against defects in materials and workmanship from date of purchase. Warranty includes the light module.

Specifications

Operation method

- Control Panel with membrane switches

Test Slide Advance/Reverse

- Via control panel

Natural Line-of-Sight Testing

- No

Comfort Height Adjustment

- Yes

Electronic Pointing System

- No

New Fluorescent Light Source

- Yes

Portability

- Integrated Handle

Soft Foam Headrest

- Yes

Head Positioning Sensor

- Yes

Peripheral Vision Testing

- Yes

Night Vision Testing

- Yes

Carrying Case (soft-sided)

- Optional

Power Supply

- 110-240 VAC, 0.4-0.2 A, 50/60 Hz

Illumination

- As per ISO 8596 and ANSI Z80.21 Standards

Dimensions

- 11" x 17" x 7.75" (28 x 43 x 19 cm), W x L x H in closed position

Weight

- 13 lbs / 6 kg

Temperature

- +50°F to +104°F (+10°C to +40°C)

Humidity

- 30% to 75%

Standard Accessories

- 8 Test Slides made from high grade photographic film sealed between optical-quality glass • Training Manual • Record Forms: pad of 100 forms to record test results • 10 Lens Cleaning Wipes • 3 Fog Eliminator Cloths • Accessory Case for storing lens cleaning wipes, fog eliminator cloths, and lenses • Dust Cover for instrument storage

Optional Accessories

- Carrying Case-soft-sided, ergonomically designed with wheels, for easy transport of vision screener • Intermediate Lenses - 19, 22, 26, 32 and 40 inches (50, 57, 67, 80 and 100 cm) • Plus Lenses with values of +1.00, +1.50, +1.75 and +2.25 are available for testing children's vision

Literature & Documents

Titmus Vision Screener V4 ASNZS

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/4294998582/1033.aspx

Manual Titmus V4, Titmus V2

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/36266/1033.aspx

Data Sheet for Occupational Landolt Model

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/4294976904/1033.aspx

Data Sheet for Occupational Letters Model

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/4294976905/1033.aspx

Data Sheet Vision Screener Professional Model_Eng

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/19279/1033.aspx

Data Sheet Vision Screener Pediatric Model_Eng

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/19276/1033.aspx

Data Sheet Vision Screener Aeromedical Model_Eng

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/19260/1033.aspx

Data Sheet Vision Screener School Model_Eng

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/19282/1033.aspx

Data Sheet Vision Screener Preschool Model_Eng

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/19278/1033.aspx

Data Sheet Vision Screener Driver Education Model_Eng

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Vision_Screening/19266/1033.aspx

Training & Proper Use

Download User Instruction Manual for the Titmus V4 and the Titmus V2 Vision Screener.

Additional Information

Email: info@titmus.com

Tel: 1-800-446-1802 (US) or +1-804-452-5200 (Intl)

Technical Support

France : info-france.hsp@honeywell.com

Export : info-export.hsp@honeywell.com

Belgium, Luxembourg, netherlands : info-benelux.hsp@honeywell.com

Germany, Austria, Switzerland : info-germany.hsp@honeywell.com

Italy : info-italia.hsp@honeywell.com

United Kingdom, Middle East & SouthernAfricaBranch : info-uk.hsp@honeywell.com

Spain, Portugal : info-spain.hsp@honeywell.com

Sweden, Denmark, Finland, Iceland, Norway : info-nordic.hsp@honeywell.com

Estonia, Latvia, Lithuania : info-nordic.hsp@honeywell.com

Russia : info-russia.hsp@honeywell.com

Hungary : info-hungary.hsp@honeywell.com

Poland : info-poland.hsp@honeywell.com

Australia/New Zealand: hshspmarketing@honeywell.com

RELATED PRODUCTS

